



All M301 orders includes:

1. Production in serum free medium if possible.
2. Purification with proteinA or proteinG or other methods depending on subclass.
3. Steril filtration of purified antibodies.
4. Measurement of protein content by UVspectroscopy with A280nm method.
5. SDS-page gel electrophorese.

Before start of production we must ensure that the cell lines are free of mycoplasmas. This is inclusive with cell lines from our lab. Using cell lines from outside we need a mycoplasm negative test from the customer or we can offer this test in our lab (Cat.No. M101).